

# PsPM: Default Parameters

May 15, 2023

Version 6.1.0

by the PsPM team<sup>1</sup>:

Dominik R Bach, Giuseppe Castegnetti, Laure Ciernik, Samuel Gerster, Saurabh Khemka, Christoph Korn,  
Samuel Maxwell, Tobias Moser, Philipp C Paulus, Ivan Rojkov, Matthias Staib, Eshref Yozdemir, Teddy  
Zhao and collaborators

---

<sup>1</sup>If you have comments on or error corrections to this documentation, please send them to the PsPM team or post them on: [bach-lab.org/pspm](https://bach-lab.org/pspm)

# Contents

<b>1</b>	<b>Introduction</b>	<b>4</b>
<b>2</b>	<b>Index</b>	<b>4</b>
2.1	Field name . . . . .	4
2.2	Data type . . . . .	4
2.3	Unit . . . . .	4
2.4	Default value . . . . .	4
2.5	Acceptable values . . . . .	4
<b>3</b>	<b>References for Values</b>	<b>5</b>
3.1	Abbreviations . . . . .	5
3.2	Action descriptors . . . . .	5
3.3	Others . . . . .	6
<b>4</b>	<b>Function-specific Default Values</b>	<b>7</b>
4.1	Blink Saccade Filter . . . . .	7
4.2	Compute Visual Angle . . . . .	7
4.3	Con1 . . . . .	7
4.4	Con2 . . . . .	7
4.5	Convert Area To Diameter . . . . .	7
4.6	Convert Au To Unit . . . . .	8
4.7	Convert ECG To Heartbeat . . . . .	8
4.8	Convert ECG To Heartbeat Amri . . . . .	8
4.9	Convert Gaze Distance . . . . .	9
4.10	Convert Heartbeat To HP . . . . .	9
4.11	Convert Pixel To Unit . . . . .	9
4.12	Convert PPG To Heartbeat . . . . .	9
4.13	Convert Visual Angle To SPS . . . . .	9
4.14	Data Editor . . . . .	10
4.15	DCM . . . . .	10
4.16	DCM Inversed . . . . .	10
4.17	Down . . . . .	11
4.18	ECG Editor . . . . .	11
4.19	EMG PP . . . . .	12
4.20	Exp . . . . .	12
4.21	Extract Segments . . . . .	12
4.22	Find Sounds . . . . .	13
4.23	Find Valid Fixations . . . . .	13
4.24	Gaze Preprocessing . . . . .	14
4.25	Get Marker Information . . . . .	14
4.26	Get RF . . . . .	14
4.27	GLM . . . . .	15
4.28	Import . . . . .	15
4.29	Interpolate . . . . .	16
4.30	Load1 . . . . .	16

4.31 Merge . . . . .	16
4.32 PFM . . . . .	16
4.33 PP . . . . .	17
4.34 Process illuminance . . . . .	17
4.35 Pupil Correct Eyelink . . . . .	17
4.36 Pupil Preprocessing . . . . .	18
4.37 Remove Epochs . . . . .	18
4.38 Resp Preprocessing . . . . .	18
4.39 SCR Preprocessing . . . . .	18
4.40 Segment Mean . . . . .	19
4.41 SF . . . . .	19
4.42 SF DCM . . . . .	20
4.43 SF MP . . . . .	20
4.44 Split Sessions . . . . .	20
4.45 Trim . . . . .	21
4.46 Write Channel . . . . .	21

# 1 Introduction

Option fields are the fields of the options that are used by many functions in PsPM. Some of them have strict specification requirements and should be set carefully if non-standard values are used. Please check this chapter carefully about how to set the fields.

## 2 Index

### 2.1 Field name

- A word written in bold means the name of a variable.
  - Example: **channel**

### 2.2 Data type

There are five basic kinds of data types used for the values of the fields in PsPM, which are namely *cell*, *character* (abbreviated as *char*), *double*, *logical*, and *struct*. Among *double*, there are some values required to be *integers* for their actual meanings. There are also some values required to be more than  $1 \times 1$  size, and there are therefore denoted as *double (vector)* or *double (matrix)*. If unspecified, *double (vector)* denotes a matrix of size  $1 \times n$ . If such vectors / matrixes are additionally required to be *integers*, they will be specified as *integer (vector)* or *integer (matrix)*.

### 2.3 Unit

Content in this column demonstrates the unit of default and acceptable values of this field. This column will not be available for the tables when there are no variables with physical meaning and units.

### 2.4 Default value

Default values are used if users have not customised the fields.

- A word written in typewriter format means a string that is used in the code, typically the value of a variable.
  - Example: add
- A number written in typewriter denotes a number used in the code, typically the value of a variable.
  - Example: 1

### 2.5 Acceptable values

Apart from the notations described above for the default values, two additional terms, Any and Subset, are also used:

Any      *Any* means any value that meets the requirement of data type can be used here.

Subset      *Subset* means any value that is a subset of the default value can be used here.

### 3 References for Values

#### 3.1 Abbreviations

ecg	Electrocardiogram
emg	Electromyography
dcm	Dynamic Causal Modelling
teo	Teager Energy Operator

#### 3.2 Action descriptors

channel\_action

- none
- add
- replace

channel\_output

- all
- corrected

eyes

- combined
- left
- right

nan\_output

- none
- screen
- file\_output

statstype

- param
- cond
- recon      Reconstruction

timeunit

- seconds
- samples
- markers

combined    The combined left and right eyes.

replace     The action of replacing the current content by the results for the channel of interest.

overwrite

1            The results will overwrite the current data.

0            The results will be dropped and not overwrite the current data.

### **3.3 Others**

NS           The parameter does not have a default value.

## 4 Function-specific Default Values

### 4.1 Blink Saccade Filter

**Table 1** Default values for pspm\_blink\_saccade\_filt

Field name	Data type	Default value	Acceptable values
channel	cell / char / integer	0	Any
channel_action	char	add	replace

### 4.2 Compute Visual Angle

**Table 2** Default values for pspm\_compute\_visual\_angle

Field name	Data type	Default value	Acceptable values
channel_action	char	add	replace

### 4.3 Con1

**Table 3** Default values for pspm\_con1

Field name	Data type	Default value	Acceptable values
zscored	logical	0	1

### 4.4 Con2

**Table 4** Default values for pspm\_con2

Field name	Data type	Default value	Acceptable values
overwrite	double	0	1, 2

### 4.5 Convert Area To Diameter

**Table 5** Default values for pspm\_convert\_area2diameter

Field name	Data type	Default value	Acceptable values
channel_action	char	add	replace

## 4.6 Convert Au To Unit

**Table 6** Default values for pspm\_convert\_au2unit

Field name	Data type	Default value	Acceptable values
channel_action	char	add	replace

## 4.7 Convert ECG To Heartbeat

**Table 7** Default values for pspm\_convert\_ecg2hb

Field name	Data type	Unit	Default value	Acceptable values
channel_action	char	/	replace	add
debugmode	logical	/	0	1
maxHR	double	bpm	200	> 20
minHR	double	bpm	20	< 200
outfact	double	/	2	Any
semi	logical	/	0	1
twthresh	double	second	0.36	Any

## 4.8 Convert ECG To Heartbeat Amri

**Table 8** Default values for pspm\_convert\_ecg2hb\_amri

Field name	Data type	Unit	Default value	Acceptable values
channel	cell / char / integer	/	ecg	Any
channel_action	char	/	add	replace
ecg_bandpass	double (vector)	bpm	[0.5, 40]	[m, n]: m>0, n>0, n>m
hrrange	double (vector)	bpm	[20, 200]	[m, n]: m>0, n>0, n>m
min_cross_corr	double	/	0.5	Any
min_relative_amplitude	double	/	0.4	Any
signal_to_use	char	/	auto	ecg, teo
teo_bandpass	double (vector)	Hz	[8, 40]	> 0
teo_order	integer	/	1	Any



## 4.9 Convert Gaze Distance

**Table 9** Default values for pspm\_convert\_gaze\_distance

Field name	Data type	Default value	Acceptable values
channel_action	char	add	replace

## 4.10 Convert Heartbeat To HP

**Table 10** Default values for pspm\_convert\_hb2hp

Field name	Data type	Default value	Acceptable values
channel_action	char	replace	add
limit.lower	double	0.2	> 0
limit.upper	double	2	> 0

## 4.11 Convert Pixel To Unit

**Table 11** Default values for pspm\_convert\_pixel2unit

Field name	Data type	Default value	Acceptable values
channel_action	char	add	replace

## 4.12 Convert PPG To Heartbeat

**Table 12** Default values for pspm\_convert\_ppg2hb

Field name	Data type	Default value	Acceptable values
channel_action	char	replace	add
diagnostics	logical	0	1

## 4.13 Convert Visual Angle To SPS

**Table 13** Default values for pspm\_convert\_visangle2sps

Field name	Data type	Default value	Acceptable values
channel	cell / char / integer	1	Any
channel_action	char	add	replace
eye	char	settings.lateral.char.b	settings.lateral.char.l, settings.lateral.char.r

## 4.14 Data Editor

**Table 14** Default values for pspm\_data\_editor

Field name	Data type	Default value	Acceptable values
epoch_file	char	NS	file must be a struct with an epoch field
output_file	char	NS	a file the changed data is saved to
overwrite	double	0	1, 2

## 4.15 DCM

**Table 15** Default values for pspm\_dcm

Field name	Data type	Default value	Acceptable values
aSCR_sigma_offset	double	0.1	> 0
method	char	dcm	?
sclpost	double	5	> 0
sclpre	double	2	> 0
sfreq	double	0.5	> 0
sfpost	double	5	> 0
sfpre	double	2	> 0
crfupdate	logical	0	1
depth	integer	2	Any
dispsmallwin	logical	0	1
dispwin	logical	1	0
eventnames	cell	{}	Any
getrf	logical	0	1
indr	logical	0	1
nosave	logical	0	1
overwrite	double	1	0, 2
rf	logical	0	1
trlnames	cell / char	{}	Any

## 4.16 DCM Inversed

**Table 16** Default values for pspm\_dcm\_inv

Field name	Data type	Default value	Acceptable values
aSCR	double	0	?
aSCR_sigma_offset	double	0.1	> 0
eSCR	logical	0	1
sclpost	double	5	> 0
sclpre	double	2	> 0
sffreq	double	0.5	> 0
sfpost	double	5	> 0
sfpre	double	2	> 0
crfupdate	logical	0	1
depth	integer	2	Any
dispsmallwin	logical	0	1
dispwin	logical	1	0
eventnames	cell	{}	Any
getrf	logical	0	1
indrf	logical	0	1
nosave	logical	0	1
overwrite	double	1	0, 2
rf	logical	0	1
trlnames	cell / char	{}	Any

#### 4.17 Down

**Table 17** Default values for pspm\_down

Field name	Data type	Default value	Acceptable values
overwrite	double	0	1, 2

#### 4.18 ECG Editor

**Table 18** Default values for pspm\_ecg\_editor

Field name	Data type	Default value	Acceptable values
artefact	char	[]	Any
channel	cell / char / integer	1	Any
factor	double	1	> 0
semi	logical	0	1

#### 4.19 EMG PP

**Table 19** Default values for pspm\_emg\_pp

Field name	Data type	Default value	Acceptable values
channel	cell / char / integer	emg	Any
channel_action	char	replace	add
mains_freq	double	50	> 0

#### 4.20 Exp

**Table 20** Default values for pspm\_exp

Field name	Data type	Default value	Acceptable values
delim	char	\t	Any
exclude_missing	logical	0	1
statstype	char	param	cond, recon
target	char	screen	Any

#### 4.21 Extract Segments

**Table 21** Default values for pspm\_extract\_segments

Field name	Data type	Default value	Acceptable values
marker_chan	cell / char / integer	marker	Any
length	double	10	$\geq 0$
nan_output	char	none	screen, file output
norm	logical	0	1
outputfile	char	‘ ’	Any
overwrite	double	0	1, 2
plot	logical	0	1
timeunit	char	seconds	samples, markers

## 4.22 Find Sounds

**Table 22** Default values for pspm\_find\_sounds

Field name	Data type	Default value	Acceptable values
channel_action	char	none	add, replace
channel_output	char	all	corrected
diagnostics	logical	1	0
expectedSoundCount	integer	0	$\geq 0$
maxdelay	double	3	$\geq 0$
mindelay	double	0	$\geq 0$
plot	logical	0	1
resample	integer	1	$\geq 1$
roi	double (vector)	[]	[a, b]; a,b $\in\mathbb{R}$
sndchannel	integer	0	$\geq 0$
threshold	double	0.1	$\geq 0$
trigchannel	integer	0	$\geq 0$

## 4.23 Find Valid Fixations

**Table 23** Default values for pspm\_find\_valid\_fixations

Field name	Data type	Default value	Acceptable values
channel	cell / char / integer	1	Any
eyes	char	combined	left, right
fixation_point	double (vector)	[0.5, 0.5]	?
missing	logical	0	1
newfile	char	‘ ’	Any
overwrite	double	0	1, 2
plot_gaze_coords	logical	0	1
resolution	double (vector)	[1, 1]	[a, b]; a,b∈R

## 4.24 Gaze Preprocessing

**Table 24** Default values for pspm\_gaze\_pp

Field name	Data type	Default value	Acceptable values
channel	char	gaze_x_l	gaze_x_r, gaze_y_l, gaze_y_r
channel_action	char	add	replace
channel_combine	char	none	gaze_x_l, gaze_x_r, gaze_y_l, gaze_y_r
valid_sample	logical	0	1

## 4.25 Get Marker Information

**Table 25** Default values for pspm\_get\_markerinfo

Field name	Data type	Default value	Acceptable values
filename	char	‘ ’	Any
marker_chan	double	-1	Any
overwrite	double	0	1, 2

## 4.26 Get RF

**Table 26** Default values for pspm\_get\_rf

Field name	Data type	Default value	Acceptable values
aSCR_sigma_offset	double	0.1	> 0
eventnames	cell	{}	Any
nosave	logical	0	1
sclpost	double	5	> 0
sclpre	double	2	> 0
sffreq	double	0.5	> 0
sfpost	double	5	> 0
sfpre	double	2	> 0
crfupdate	logical	0	1
depth	double	2	?
dispsmallwin	logical	0	1
dispwin	logical	1	0
getrf	logical	0	1
indr	logical	0	1
overwrite	double	1	0, 2
rf	logical	0	1
trlnames	char / cell	{}	Any

## 4.27 GLM

**Table 27** Default values for `pspm_glm`

Field name	Data type	Default value	Acceptable values
marker_chan_num	cell / char / integer	marker	Any
bf	logical	0	1
exclude_missing	struct	NS	struct('segment_length', m, 'cutoff', n), m, n > 0
centering	logical	1	0
norm	logical	0	1
overwrite	double	0	1, 2

## 4.28 Import

**Table 28** Default values for `pspm_import`

Field name	Data type	Default value	Acceptable values
overwrite	double	0	1, 2

## 4.29 Interpolate

**Table 29** Default values for pspm\_interpolate

Field name	Data type	Default value	Acceptable values
channel	cell / char / integer	1	Any
channel_action	char	add	replace
extrapolate	logical	0	1
method	char	linear	pchip, nearest, spline, previous, next
newfile	logical	0	1
overwrite	double	0	1, 2

## 4.30 Load1

**Table 30** Default values for pspm\_load1

Field name	Data type	Default value	Acceptable values
overwrite	integer	0	1, 2
zscored	logical	0	1

## 4.31 Merge

**Table 31** Default values for pspm\_merge

Field name	Data type	Default value	Acceptable values
marker_chan_num	integer (vector)	[0, 0]	Any
overwrite	integer	0	1, 2

## 4.32 PFM

**Table 32** Default values for pspm\_pfm

Field name	Data type	Default value	Acceptable values
overwrite	integer	0	1, 2



### 4.33 PP

**Table 33** Default values for pspm\_pp

Field name	Data type	Default value	Acceptable values
overwrite	integer	0	1, 2

### 4.34 Process illuminance

**Table 34** Default values for pspm\_process\_illuminance

Field name	Data type	Default value	Acceptable values
bf	struct	struct()	Any
fn	char	<i>empty</i>	Any
overwrite	integer	0	1, 2
transfer	double (vector)	[49.79, -1.05, -0.50]	[a, b, c]: a, b, c>0
bf.constriction	struct	struct()	Any
bf.dilation	struct	struct()	Any
bf.duration	double	20	$\geq 0$
bf.offset	double	0.2	$\geq 0$

### 4.35 Pupil Correct Eyelink

**Table 35** Default values for pspm\_pupil\_correct\_eyelink

Field name	Data type	Default value	Acceptable values
C_x	double	0	Any
C_y	double	0	Any
C_z	double	0	Any
channel	char	pupil	Any
channel_action	char	add	replace
mode	char	auto	manual
S_x	double	0	Any
S_y	double	0	Any
S_z	double	0	Any
screen_size_mm	double (vector)	[43.5, 29.9]	[a, b]: a, b>0
screen_size_px	double (vector)	[1920, 1080]	[a, b]: a, b>0

## 4.36 Pupil Preprocessing

**Table 36** Default values for `pspm_pupil_pp`

Field name	Data type	Default value	Acceptable values
<code>channel</code>	char	<code>pupil</code>	<code>pupil_l</code> , <code>pupil_r</code>
<code>channel_combine</code>	char	<code>none</code>	<code>pupil_l</code> , <code>pupil_r</code>
<code>plot_data</code>	logical	<code>0</code>	<code>1</code>
<code>segments</code>	cell	<code>{}</code>	Any

## 4.37 Remove Epochs

**Table 37** Default values for `pspm_remove_epochs`

Field name	Data type	Default value	Acceptable values
<code>channel_action</code>	char	<code>add</code>	<code>replace</code>

## 4.38 Resp Preprocessing

**Table 38** Default values for `pspm_resp_pp`

Field name	Data type	Default value	Acceptable values
<code>channel_action</code>	char	<code>add</code>	<code>replace</code>
<code>datatype</code>	cell	<code>{rp, ra, rfr, rs, all}</code>	Subset
<code>diagnostics</code>	logical	<code>0</code>	<code>1</code>
<code>plot</code>	logical	<code>0</code>	<code>1</code>
<code>systemtype</code>	char	<code>bellows</code>	<code>cushion</code>

## 4.39 SCR Preprocessing

**Table 39** Default values for `pspm_scr_pp`

Field name	Data type	Default value	Acceptable values
change_data	logical	1	0
channel_action	char	add	replace, withdraw
clipping_n_window	integer	10000	Any
clipping_step_size	integer	2	Any
clipping_threshold	double	0.1	Any
data_island_threshold	double	0	$\geq 0$
deflection_threshold	double	0.1	Any
expand_epochs	double	0.5	$\geq 0$
max	double	60	$> 0$
min	double	0.05	$> 0$
missing_epochs_filename	char	missing_epochs_filename	Any
slope	double	10	Any

#### 4.40 Segment Mean

**Table 40** Default values for pspm\_segment\_mean

Field name	Data type	Default value	Acceptable values
adjust_method	char	none	downsample, interpolate
newfile	char	<i>empty</i>	Any
overwrite	integer	0	1, 2
plot	logical	0	1

#### 4.41 SF

**Table 41** Default values for pspm\_sf

Field name	Data type	Default value	Acceptable values
dispsmallwin	logical	0	1
dispwin	logical	1	0
fresp	double	0.5	$\geq 0$
marker_chan_num	char / integer	marker	Any
overwrite	integer	1	0, 2
theta	double (vector)	[0.92, 3.92, 2.16, 1.53, 1.64]	Any
threshold	double	0.1	$> 0$

#### 4.42 SF DCM

**Table 42** Default values for pspm\_sf\_dcm

Field name	Data type	Default value	Acceptable values
dispwin	logical	1	0
dispsmallwin	logical	0	1
fresp	double	0.5	$> 0$
theta	double (vector)	[0.92, 3.92, 2.16, 1.53, 1.64]	Any
threshold	double	0.1	$> 0$

#### 4.43 SF MP

**Table 43** Default values for pspm\_sf\_mp

Field name	Data type	Default value	Acceptable values
diagnostics	logical	0	1
dispwin	logical	0	1
fresp	double	0.5	$> 0$
theta	double (vector)	[0.92, 3.92, 2.16, 1.53, 1.64]	Any
threshold	double	0.1	$> 0$

#### 4.44 Split Sessions

**Table 44** Default values for pspm\_split\_sessions

Field name	Data type	Default value	Acceptable values
max_sn	double	settings.split.max_sn	> 0
min_break_ratio	double	settings.split.min_break_ratio	> 0
missing	char	<i>empty</i>	Any
overwrite	integer	0	1, 2
prefix	double	0	$\leq 0$
randomITI	logical	0	1
splitpoints	double (vector)	[]	Any
suffix	double	0	$\geq 0$
verbose	logical	1	0

#### 4.45 Trim

**Table 45** Default values for pspm\_trim

Field name	Data type	Default value	Acceptable values
drop_offset_markers	integer	0	Any
marker_chan_num	integer	0	Any
overwrite	integer	0	1, 2

#### 4.46 Write Channel

**Table 46** Default values for pspm\_write\_channel

Field name	Data type	Default value	Acceptable values
channel	integer / char / cell	0	Any
delete	char	last	first, all
msg	char / struct	<i>empty</i>	Any
prefix	char	Generic undocumented operation ::	Any